

pecaplasty[®]

Percutaneous Bunion Correction


Sales Brochure




Translation of the metatarsal head



Bunion deformity




Initial placement of pecaplasty



Translation of metatarsal head by turning dial

Accurate & adjustable K-wire placement




Adjustable entry point: ●

The insertion point of the screw in the metatarsal must be as proximal as possible.

Fixed aiming point: ●

Sliding of the arc around its center of rotation.



Fixed position in the medial side of the head.



Fixed distance between the aiming point and the paddle of the PECA guide.

References

PECA Targeting Guide

Reference	Designation	Image
XMS01038-1	PECA targeting guide	
XMS01038-2	PECA targeting guide wire handle	
SKW06001	PECA targeting guide wire	
XMS01038-3	PECA Ø3mm sleeve	
XMS01038-4	PECA Ø4mm sleeve	

Distributed by:

Novastep Inc.
30 Ramland Road, Suite 200
Orangeburg, NY 10962
(877) 287-0795
cs@novastep-us.com

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a surgeon.
Rx only.

Novastep Inc., or its affiliates, own, use or have applied for the following trademarks or service marks: Pecaplasty, Novastep. All other trademarks are trademarks of their respective owners or holders.

© 2021 Novastep